

## FEATURES

Large 1/2" L.E.D. displays the time and aids setting alarms.



# MedReady™

Automated Medication Dispenser  
Available with cellular  
compliance monitoring

- Medication regimen can dose up to 4 times per day:
  - ◆ "Time for medicine" alarm beeps 5 seconds on, 1 second off continuously until the medication door is slid open and the "off" position is activated.
  - ◆ Note: Disabling alarm requires the medication door to be slid open.
  - ◆ Programmable alarm sounds for 3, 10, 20, 30, ...250 minutes, as user selects, and can chirp each minute thereafter until the dose is taken or the next one is available.
  - ◆ Powerful, 90dB, 400Hz. audible alarm is designed for maximum effectiveness.
  - ◆ Audible and visual alarms can be independently turned on or off.
- Only the most recent medication dose is available. Pill tray rotates to the next dose when the next audible alarm sounds. If a patient misses a dose, the pill dispenser denies access to the missed dose.
- MedReady holds 28 doses of dry medication such as pills or capsules.
  - ◆ 4 Doses a day for 7 days
  - ◆ 3 Doses a day for 9 days
  - ◆ 2 Doses a day for 14 days
  - ◆ 1 Dose a day for 28 days
- Alarm in pill dispenser is as easy to set as a digital clock.
- Medication is under lock and key, and is tamper resistant.
- Pill compartments are "ramped" for easy removal and can fit up to 9 regular **m&m**s.
- Compartment divisions keep medications in place even if the unit is turned upside-down.
- MedReady should remain plugged into electric power, although there is a 48 hour rechargeable battery back-up for travel or power failures. MedReady should not be used in wet or very damp environments.
- "OFF" feature for long term disabling of the dispenser.
- "Early dosage" feature allows removal of medication prior to scheduled times.

